

5,000 g x 1 g


DIGITAL SACLE

Model : GM-5001

ISO-9001, CE, IEC1010

www.yalab.com.tw

02-2389-0101



LITRON

www.YaLAB.com.tw

The Art of Measurement

5,000 g x 1 g

DIGITAL SACLE

Model : GM-5001

www.YaLAB.com.tw 02-2389-0101

FEATURE

- * Microprocessor circuit with high accuracy & better performance.
- * High resolution & wide measuring capacity, 5,000 g x 1 g.
- * Large LCD display, 18 mm (0.7") digit size, easy read out.
- * Built-in gram & oz display unit, select by the internal slide switch.
- * GM-501 built-in self calibration system, high precision.
- * Heavy duty ABS housing plastic case.
- * LCD, 18 mm (0.7") digit size, 5 digits with annunciator.

SPECIFICATION

Measuring Capacity	5,000 g	
Resolution	1 g	
Min. Display Weight	3 g	
Unit Select	g or oz, select by internal slide switch.	
Accuracy (23± 5 °C)	± (0.3 % + 1 d) <i>* After the moment when the self-calibration be executed.</i>	
Sampling time	Approx. 0.8 second.	
Tare Control	Full capacity.	
Transducer	Load cell.	
Circuit	Microprocessor circuit.	
Operating Temperature	0 °C to 50 °C (32 °F to 122 °F).	
Operating Humidity	Max. 80% RH.	
Size	Cabinet	250 x 188 x 70 mm
	Platform	185 x 185 mm.
Power Supply	6 x 1.5V AA (UM-3) BATTERY, or DC 9V adapter.	
Power Consumption	Approx. DC 27 mA.	
Accessory Included	Operation Manual.....1 PC.	

* Appearance and specifications listed in this brochure are subject to change without notice.

GM-5001+1109

www.YaLAB.com.tw 02-2389-0101